On-Call® MultiPro

Kit de test RAC Notice d'emballage

REF	G136-132	#	OGS-224	Français
-----	----------	---	---------	----------

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




BIBLIOGRAPHIE

- Shaw A.B., Risdon P., Lewis-Jackson J. Protein creatinine index and Albustix in assessment of proteinuria. British Medical Journal. 287(1983):929-932.
- Rodriguez-Thompson D., Lieberman E.S. Use of a random urinary protein-to-creatinine ratio for the diagnosis of significant proteinuria during pregnancy. Am J Obstet Gynecol. 185(2001):808-811.
- Mogensen C.E.et al. Prevention of diabetic renal disease with special reference to Microalbuminuria. The Lancet. 346(1995):1080-1084.
- Row D.J.F., Dawnay A., Wattst G.F. Microalbuminuria in diabetes mellitus: review and recommendations for the measurement of albumin in urine. Ann Clin Biochem. 27(1990):297-312.
- Harmonen A., Vuorinen P., Jokela H. Turbidimetric measurement of Microalbuminuria. Clinica Chimica Acta. 166(1987):85-89.

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For testing CRP in human whole blood, serum and plasma.
For professional *in vitro* diagnostic use only.

INTENDED USE

The *On Cal® MultiPro* CRP Test Kit is to be used with the *On Cal® MultiPro* Analyzer for the quantitative determination of C-reactive protein (CRP) in human whole blood, serum and plasma. The *On Cal® MultiPro* CRP Test Kit is for professional *in vitro* diagnostic use only.

C-reactive protein (CRP) is mainly used as a non-specific inflammatory marker. A full-range CRP assay measures both conventional CRP and high-sensitivity CRP (hs-CRP) levels. Conventional CRP levels can be used in the evaluation of infection, tissue damage and inflammatory diseases, providing information for the diagnosis, treatment and monitoring of such conditions, while hs-CRP is a sensitive indicator for distinguishing low-level inflammatory states and is most commonly used as a cardiovascular aid in disease risk identification. Used in conjunction with the traditional clinical diagnosis of acute coronary syndrome, it is considered to be an early warning indicator of coronary artery disease or a recurrence of acute coronary syndrome.

SUMMARY


C-reactive protein is an acute phase protein synthesized by hepatocytes. It exists in low concentrations in healthy individuals; however, in disease states, elevated levels are associated with invasive bacterial infection, inflammation or tissue destruction. CRP levels rise rapidly within the first 6 to 8 hours of onset.¹ Likewise, they can rapidly decline, with an elimination half-life estimated at 4 to 9 hours. CRP has been identified as a useful marker for monitoring flare-ups of chronic inflammatory conditions, such as rheumatoid arthritis, inflammatory bowel disease, and a number of vasculitic syndromes. Additionally, measurement of CRP levels can be useful in the stratification of patients at intermediate risk for a cardiovascular event.²

PRINCIPLE AND REFERENCE VALUES

Anti-CRP antibodies coated on the surface of the latex particles react with the antigen (CRP) in the sample to form an antigen/antibody complex. This reaction causes changes in the turbidity of the solution and alters light absorbance, which is then measured by the *On Cal® MultiPro* Analyzer. The latex helps amplify the reactive signal.

The agglutination reaction described above causes increased scattering of light, measured as an increase in absorbance. This change in absorbance is inversely proportional to the concentration of CRP in the sample. The higher the concentration of CRP in the sample, the greater the change of absorbance during the reaction.

The calibration equation is stored in the QR code on the cartridge, which is automatically scanned by the analyzer at the beginning of each test. For whole blood samples, the test results are adjusted automatically according to hematocrit (HCT) levels.

REAGENTS

The concentrations provided below may vary within manufacturing tolerances.

0.075% w/v CRP-specific monoclonal antibody-latex; 1% w/v polyethylene glycol; 0.05% w/v saponin; 98.875% w/w nonreactive ingredients.

PRECAUTIONS

- Use only *On Cal® MultiPro* CRP Test Kits with *On Cal® MultiPro* Analyzer.
- Keep test cartridges in their sealed foil package until use.
- Do not use the test kit past the expiration date printed on the box.
- Do not touch the optical window of the test cartridges.
- Discard any contaminated or damaged test cartridges.
- Do not reuse test cartridges.
- Discard if the desiccant bag is missing or open.
- Discard if the foil package is open for more than 120 minutes.
- Specimens should be considered potentially hazardous. Handle in the same manner as an infectious agent.
- Used test cartridges should be discarded according to local regulations after testing.

STORAGE AND STABILITY

- All unopened test cartridges will remain stable through the expiration date printed on the label if stored between 36-86 °F (2-30 °C).
- Allow foil pouch containing test cartridge to reach 59-89.6 °F (15-32 °C) before opening the pouch.
- Test cartridge must be used within 120 minutes of being removed from foil pouch.
- DO NOT expose test cartridge to strong light during storage or testing.
- DO NOT freeze.

SPECIMEN COLLECTION AND PREPARATION

- Accepted sample types include fresh capillary, venous whole blood, plasma and serum.
- Compatible anticoagulants include EDTA and heparin.
- The hematocrit range of whole blood specimens is 20%-70%. The analyzer will display 'Hb Low' or 'Hb High' if the level of hematocrit is outside this range.
- After the provided Capillary Device is filled with whole blood (5µL capacity), analysis must begin within 1 minute if no anticoagulant has been added to the sample. If anticoagulant has been added to the sample, analysis must begin within 3 minutes.
- EDTA and heparin preserved whole blood should be tested within 24 hours.
- Plasma and serum may be stored at 2-8 °C (36-46 °F) for one week, or below -20 °C (≤-4 °F) for six months.
- Allow samples to reach room temperature (59-89.6 °F, 15-32 °C) and mix thoroughly before use.

MATERIALS

Materials Provided

- Test Cartridges
- Capillary Devices
- Package Insert

Materials Required But Not Provided

- Safety Lancets
- Analyzer
- Gauze for Puncture Site
- Gloves
- Alcohol Swab
- Lint-free Tissue

DIRECTIONS FOR USE

If testing materials are stored at 2-8 °C, please place test cartridge, specimen, controls, and other testing components at room temperature (15-32 °C) for at least 1 hour before testing to ensure the materials reach the required operating temperature (15-32 °C). Refer to *On Cal® MultiPro* System User's Manual for detailed instructions.

- Switch on the Analyzer.
- Remove the test cartridge from the foil pouch.
 - Note:** Discard the test cartridge if it is damaged, the desiccant is missing, or if loose desiccant particles are found inside the foil pouch.
 - After opening the foil pouch, the test cartridge must be used within 120 minutes. Do not use a test cartridge if it falls on the floor after being removed from the foil pouch.
- Remove the Capillary Device from the test cartridge. Fill glass capillary with whole blood from finger or blood collection tube, serum or plasma from collection tube. Keep the glass capillary in contact with the sample for another 2 seconds after the sample has filled the capillary device to ensure the capillary is completely full (but not spilling over).
 - Note:** DO NOT submerge the Capillary Device in the collection tube. Use a lint-free tissue to carefully wipe the outside of the glass capillary if excess sample is on the Capillary Device. Avoid touching the tip of capillary tube with the lint-free tissue. If the tissue touches the tip, and there appears to be less blood left in the tube, discard the Capillary Device and repeat the procedure using a new Capillary Device.
 - Inspect the glass capillary for the presence of bubbles. If bubbles are present, discard the Capillary Device and repeat the procedure using a new Capillary Device.
- Insert capillary device into the test cartridge vertically and slowly. Avoid scratching inner wall of cartridge with the capillary tip.
 - Note:** DO NOT touch the circular optical window of the test cartridge.
- Click the Test Settings button on the Home Screen and select the sample type on the screen. Then click Test button on the Home Screen and enter Patient ID/Patient Name/Sample ID.
- Open the analyzer cartridge compartment cover. Insert the test cartridge into the cartridge compartment with the label facing to the left. The analyzer will read the calibration equation automatically from the QR code printed on the cartridge.
 - Note:** Push the cartridge completely to the bottom of the cartridge compartment.
- Close the cover and the assay will begin automatically. The results will be displayed on the screen within 4 minutes.
 - Note:** Don't open the cover during the testing process, otherwise, the test will be canceled by the analyzer. The last step in test process is to close the cover. If the cover is closed before the Test button is clicked, the analyzer will not start the test. If this occurs, open and then close the cover to begin testing.
- Open the cartridge compartment cover. Remove the used test cartridge and close the cover.
 - Note:** Gently push the tab on the cartridge to the right when removing the cartridge.

INTERPRETATION OF RESULTS

The assay has a reportable range of: 0.5mg/L to 200mg/L CRP (0.05mg/dL-20mg/dL) in blood mode and 0.5mg/L to 120mg/L in serum and plasma mode. The test is linear throughout this range. Results lower than 0.1mg/L will be displayed as <0.1mg/L. Results higher than 200mg/L will be displayed as >200mg/L in blood mode. Results higher than 120mg/L will be displayed as >120mg/L in serum and plasma mode.

All laboratory tests are subject to random error. If you observe unexpected or questionable results, take the following steps:

- Ensure that the test cartridge is not expired.
- Compare results to controls with known levels.
- Repeat the test using a new test cartridge.
- If the problem persists, discontinue using the test cartridge. Contact your local distributor.

QUALITY CONTROL

Before test kits are released by the manufacturer, each lot undergoes a thorough analysis and characterization. Values of calibration parameters are traced to reference ERM-DA474/IFCC. To ensure quality of both testing procedures and patient results for CRP, the *On Cal® MultiPro* system performs optical, electronic, mechanical, and reagent systems checks, including calibration verification, during the course of each assay. If an assay or system error occurs during any individual measurement, the system automatically reports an error message, preventing the reporting of erroneous patient results.

For best results, performance of the CRP test kit should be confirmed by testing known specimens/controls whenever a new lot of test kits are first opened. Please check the *On Cal® MultiPro* System User's Manual for Quality Control procedures.

Each laboratory should establish its own goals for adequate standards of performance. Contact your local distributor for information on specific controls for this product.

PERFORMANCE CHARACTERISTICS

Linearity

In each sample style six test samples were analyzed in replicates of three. The linear regression is shown below:

Sample style	Linearity Equation	r
blood	y = 1.0341x + 0.1254	0.9999
serum	y = 0.9588x + 0.2166	0.9997
plasma	y = 1.0167x - 0.179	0.9999

Reproducibility and Precision

Within-run precision was evaluated by testing three serum specimens in replicates of ten on three lots of test kits.

Specimen	Specimen 1			Specimen 2			Specimen 3		
Precision	Lot 1	Lot 2	Lot 3	Lot 1	Lot 2	Lot 3	Lot 1	Lot 2	Lot 3
Mean (mg/L)	3.0	3.1	3.0	10.0	10.1	10.1	49.4	52.6	50.8
CV (%)	4.6%	5.6%	5.2%	4.5%	4.6%	4.1%	3.0%	2.9%	3.1%

Total precision is listed below:

Test level	Specimen 1	Specimen 2	Specimen 3
Mean (mg/L)	3.0	10.1	50.9
CV (%)	5.3	4.4	3.9

Accuracy

The *On Cal® MultiPro* Analyzer and *On Cal® MultiPro* CRP Test Kit were used by trained technicians to test 110 EDTA preserved venous blood specimens, 103 fresh capillary specimens, 105 serum specimens and 109 plasma specimens. The measured values were compared to an automatic biochemistry analyzer. Linear regression results are as follows:

Specimen	Slope	Intercept	R	N
Venous blood	0.9515	-0.1691	0.9982	110
Fresh capillary	1.0021	-0.0507	0.9986	103
Serum	0.9151	0.2776	0.9981	105
Plasma	0.9973	-0.4396	0.9981	109

LIMITATIONS

- The test results are for clinical reference only, and cannot be used alone as a basis for confirming or excluding cases.
- Do not use hemolyzed serum and plasma samples. Obviously turbid samples should not be used for analysis and need to be re-sampled or centrifuged.
- Be sure the sample type of the specimen is the same as the sample type selected in the analyzer test settings. Otherwise the analyzer will display 'Hb Low' or 'Hb High' instead of the test result.
- Sample dilution
Do not dilute whole blood samples.
Serum and plasma samples may be diluted with a 0.9% NaCl solution, and the dilution ratio should be less than 6:1.

The following substances will not affect test results in the concentrations listed below:

Substance	Concentration
Bilirubin	342 µmol/l
Cholesterol	25.75 mmol/L
Triglycerides	9.25 mmol/L
Rheumatoid factor	1200 IU
Leukocytes	41×10 ⁹ /L

BIBLIOGRAPHY

- Clyne B, Oshaker J S. The C-reactive protein [J]. Journal of Emergency Medicine, 1999,17(6):1019-1025
- Drs. Stephen K, Emanuele D, Lisa P. C-Reactive Protein, Fibrinogen, and Cardiovascular Disease Prediction[J]. The New England Journal of Medicine, 2012,367:1310-20

Index of Symbols

	Consult instructions for use or consult electronic instructions for use		Contents sufficient for <math>n> tests		Manufacturer
	Temperature limit		Batch code		Catalogue number
	Do not reuse	# icon"/>	Model number		Use-by date
	Do not use if package is damaged and consult instructions for use		Authorized representative in the European Community / European Union		In vitro diagnostic medical device

ACON®
ACON Laboratories, Inc.
5850 Oberlin Drive, #340
San Diego, CA 92121, USA
www.aconlabs.com

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MDSS GmbH