REF	G136-142	#	OGS-223	Français
-----	----------	---	---------	----------

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REF G126-182

English

For validating ACR testing.
For professional *in vitro* diagnostic use only.

INTENDED USE

The *On Cal® MultiPro* ACR Controls are intended for use with the *On Cal® MultiPro* Analyzer and *On Cal® MultiPro* ACR Test Kit as a quality control check. They provide confirmation that the test kit and analyzer are working properly and ensure that the test is being performed correctly. The *On Cal® MultiPro* ACR Controls are for professional *in vitro* diagnostic use only.

PRECAUTIONS

- For *in vitro* diagnostic use only. Use before the expiration date printed on the package label.
- Make sure the Controls and all test materials reach an operating temperature of 15-32°C (59-89.6°F) prior to testing. Temperatures outside that range may lead to incorrect results.
- All materials should be considered potentially hazardous and handled in the same manner as chemical agents.
- This product is not intended for use as a calibrator.
- Use the Controls before urine testing as an objective method to assess the techniques or practices in use.
- Use the *On Cal® MultiPro* ACR Controls with the *On Cal® MultiPro* Analyzer and *On Cal® MultiPro* ACR Test Kit.
- Discard the Controls if the reconstitution solution package seal is damaged or the lyophilized material is missing from the bottles.
- Controls and the Control ID Cartridge should have the same batch code (lot number).
- Gently push tab on Control ID cartridge to the right when removing the cartridge from the analyzer cartridge compartment.
- The used Test Kit and Controls should be discarded according to local regulations after testing.

COMPONENTS

The *On Cal® MultiPro* ACR Controls are available in two levels. The main components of the control products include human albumin and creatinine, Tris buffer >99% w/w and preservative proclin300 ≤ 0.002%.

STORAGE AND STABILITY

- All unopened Controls are stable through the expiration date printed on the box if stored between 36-46 °F (2-8 °C). If stored at room temperature (up to 86 °F/30 °C), unopened Controls will remain stable for 3 months unless the expiration date printed on the label has been reached.
- Controls should be reconstituted within 30 minutes after opening the bottles.
- Tighten the cap on each control bottle immediately after reconstitution or use. Store bottles in the control box between 36-46 °F (2-8 °C).
- Reconstituted Controls will remain stable for 3 months.
- DO NOT use after expiration date.
- DO NOT expose Controls to strong light during storage or testing.
- DO NOT freeze after reconstitution.

MATERIALS REQUIRED

- Control Level 1
- Control Level 2
- Control ID Cartridge
- Package Insert
- Reconstitution Solution
- White Dropper Cap
- Black Dropper Cap

MATERIALS REQUIRED BUT NOT PROVIDED

- On Cal® MultiPro* Analyzer
- On Cal® MultiPro* ACR Test Kit

INSTRUCTIONS FOR USE

If testing materials are stored at 2-8 °C, please place test cartridge, specimen, controls, and other testing components at room temperature (15-32 °C) for at least 1 hour before testing to ensure the materials reach the required operating temperature (15-32 °C). Refer to *On Cal® MultiPro* System User's Manual for detailed instructions.

Reconstitution: The following reconstitution procedures are recommended to minimize variation resulting from different methods in different laboratories.

- Gently tap the bottom of the control bottle on the counter to collect as much lyophilized material as possible on the bottom of the bottle.
- Carefully remove the cap from the control bottle.
- Open the aluminum foil bag containing the reconstitution solution. Take out the sealed plastic tube of reconstitution solution, and gently tap the tube 3-5 times to ensure that the reconstitution solution is at the bottom of the tube. Unscrew the front seal of the tube. Then add all the reconstitution solution to the control bottle.
Note: All the reconstitution solution should be added to the control bottle and not spilled.
- Tighten the cap of the control bottle and let stand at room temperature for 10 minutes. Then mix the control bottle by rotating the bottle until the solution is homogeneous and all lyophilized material is reconstituted.
Note: No foam should appear during the mixing process! Don't shake it!
- Replace the cap of Control Level 1 with a white dropper cap and replace the cap of Control Level 2 with a black dropper cap. Be sure to tighten both dropper caps.

Run a Control:

- Switch on the Analyzer.
- Select QC Test on the Home screen to enter QC Test mode. Open the cartridge compartment cover. Insert the Control ID Cartridge into the cartridge compartment with the label facing left. Close the compartment cover. The analyzer will automatically read the control information from the QR code printed on the Control ID Cartridge.
- Open the cartridge compartment cover and remove the Control ID Cartridge.
- Select Control Level 1 and/or Control Level 2 according to the prompts displayed on the analyzer screen. The concentration level of each Control is shown on the bottle label.
- Unwrap a Test Cartridge and remove the empty Capillary Device from the Test Cartridge.
- Unscrew the dropper cap. While applying only slight pressure to the bulb, insert the tip of dropper into the control solution (tilt bottle as necessary). Release pressure on bulb to aspirate a very small amount of control solution.
- Use the Capillary Device to absorb the control solution collected in the dropper and completely fill the Capillary Device.
- Do not touch the dropper to any other surfaces. Squeeze any excess control solution out of the dropper back into the control solution bottle. Carefully screw the dropper cap back onto control bottle.
- Insert filled Capillary Device back into the Test Cartridge vertically and slowly. Avoid scratching inner wall of cartridge with the capillary tip. Open the cartridge compartment cover. Insert the Test Cartridge into the cartridge compartment with the label facing left. Click Test button on the screen. The analyzer will read the calibration equation automatically from the QR code printed on the cartridge.
Note: Use a lint-free tissue to carefully wipe the outside of the glass capillary if excess control solution is on the capillary device. Avoid touching the tip of glass capillary tube with the lint-free tissue. If the tissue touches the tip, and there appears to be less control solution left in the tube, discard the capillary device and repeat the procedure using a new capillary device.
Inspect the glass capillary tube for the presence of bubbles. If bubbles are present, discard the capillary device and repeat the procedure using a new capillary device.
DO NOT touch the circular optical window of the test cartridge. Push the cartridge completely to the bottom of the cartridge compartment.
- Close the cover to start the test. Results will be displayed on the screen within 6.5 minutes.
Note: The last step in QC test process is to close the cover. If the cover is closed before the Test button is clicked, the analyzer will not start the test. If this occurs, open and then close the cover to begin testing.
- Open the cartridge compartment cover. Remove the used Test Cartridge and close the cartridge compartment cover.

EXPECTED RESULTS

The results should fall within the range printed on the bottle label and are specific for each lot of controls. If the results for each control fall within the specified control range, the *On Cal® MultiPro* system is working correctly and the procedures are being performed properly.

If the results do not fall within the specified range:

- Check the expiration date of the Test Cartridge and the Controls. Make sure they are within the expiration date. Discard any expired Test Cartridges and Controls.
 - Make sure the Controls were reconstituted within the past 3 months.
 - Make sure the Controls were fully mixed before they were run.
 - Confirm that you are using the *On Cal® MultiPro* ACR Controls.
 - Make sure that all test procedures were followed correctly.
 - Make sure the analyzer is not contaminated.
 - Make sure the Control was tested using the QC Test mode on the analyzer.
 - Make sure the Controls and Control ID Cartridge have the same batch code (lot number).
 - Make sure the operating temperature is 15-32°C (59-89.6°F).
- Note:** Allow the Controls and Test Cartridge to reach room temperature for 30 minutes if stored at 2-8 °C (36-46 °F).

After checking all of the above, please repeat the test with a new Test Cartridge. If the results still fall outside the specified control range, the analyzer may not be working properly. Contact your local distributor for further assistance.

INDEX OF SYMBOLS

 IVD	<i>In vitro</i> diagnostic medical device	 CTRL	Control range	 2°C	Temperature limit
 LOT	Batch code		Use-by date	 REF	Catalogue number
	Consult instructions for use or consult electronic instructions for use		Manufacturer		Do not use if package is damaged and consult instructions for use
	Health hazard	 EC REP	Authorized representative in the European Community / European Union		

ACON Laboratories, Inc.
5850 Oberlin Drive, Suite 340,
San Diego, CA 92121, USA
www.aconlabs.com

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MDSS GmbH